Cloud-Clone Insights: Pre-treatment Guidance for Special Biological Samples in ELISA Testing

Practical tips for handling complex specimens to reduce matrix interference and preserve target-protein activity

HUSTON, TX, UNITED STATES, August 25, 2026 /EINPresswire.com/ -- ELISA is a rapid, sensitive and reliable quantitative analytical technique. Appropriate sample handling is critical to successful ELISA measurements. Collection timing, manipulation procedures and storage conditions all directly shape final assay results. Previous technical content covered routine specimens such as blood, tissue homogenates, cell lysates and cell-culture supernatants. This follow-up piece shares hands-on handling suggestions for less-common, complex biological materials, helping researchers worldwide minimise experimental artefacts and acquire credible ELISA data.

1. Sputum samples

1.Weigh the viscous portion of collected sputum. Add 0.1 % dithiothreitol (DTT) solution at twice the volume of sputum to dissolve mucus. Repeatedly pipette the mixture, vortex for 15 seconds, then shake in a 37 °C water bath for 5 minutes.

2.Add PBS buffer equal to twice the original sputum volume and continue shaking for 15-20 minutes. Filter through a 150-mesh metal sieve, centrifuge at 1500 rpm for 10 minutes, and collect supernatant for testing.

2. Faecal samples

Collect relatively solid faecal material; overly watery samples are difficult to process and may compromise measurement accuracy. Obtain more than 50 mg of faecal specimen. Wash the sample three times with PBS, following a final mass-to-volume ratio of faeces to PBS of 1:9. Perform sonication (or mechanical homogenisation), then centrifuge at 5000 × g for 10 minutes and retain supernatant for detection.

3. Saliva samples

Place collected saliva inside centrifuge tubes and freeze at −70 °C for one hour. Thaw samples on ice, centrifuge at 2000 × g for 10 minutes and collect supernatant for testing.

4. Seminal plasma samples

Collect semen in sterile containers. Fresh semen normally appears thick and gelatinous. Allow complete liquefaction at room temperature or inside a 37 °C water-bath, driven by endogenous prostate fibrinolytic enzymes. Once fully liquefied, centrifuge semen at 4000 rpm for

10 minutes to separate seminal plasma for subsequent testing.

5. Skin tissue samples
Homogenise 20-50 mg skin tissue with pre-chilled PBS using a tissue homogeniser for full lysis. Centrifuge homogenate at 10 000 g for 20 minutes and collect supernatant for ELISA measurement.

6. Milk samples
Centrifuge milk at 12 000 g for 30 minutes under 4 °C conditions. Discard upper-layer fat and bottom-layer sediment, retain the middle liquid fraction for testing.

7. Cerebrospinal fluid and broncho-alveolar lavage fluid
After specimen collection, centrifuge at 1000 × g for 20 minutes. Harvest supernatant and store samples at −20 °C. Repeated freeze-thaw cycles must be avoided.

General sample-storage key notes
Sample processing aims to preserve target proteins. Since proteins are prone to denaturation and degradation, all manipulation steps should be carried out as gently as possible. Post-processing storage strongly impacts data quality:
Aliquot and seal processed samples to prevent repeated freeze-thaw events.
Storage duration limits: less than 1 week at 4 °C; maximum 1 month at −20 °C; maximum 2 months at −80 °C.
Allow specimens to slowly equilibrate to room-temperature before use; never rapidly thaw samples via heating, which causes protein inactivation and analyte breakdown.

Closing Remarks
Proper specimen handling represents the very first step toward dependable ELISA results. Complex non-routine biological samples bring extra challenges for researchers in diverse research fields ranging from respiratory studies to reproductive biology. By adopting gentle processing workflows and strictly following cold-storage specifications, scientists can effectively protect protein integrity and lower matrix-derived interference. Cloud-Clone will keep publishing practical technical insights, supporting laboratories around the globe to achieve stable and reproducible immunoassay outcomes.

About Cloud-Clone Corp.
Cloud-Clone Corp. is dedicated to the development and production of high-quality immunoassay reagents and detection solutions. With a focus on antibody engineering, multiplex assay development, and cross-platform compatibility, the company provides research tools designed to support precision medicine and advanced biomedical investigation globally. Our core products and services include the research and development of proteins, antibodies, ELISA kits, primary cells, and multiplex cytokine assay kits, as well as professional CRO services to fully meet the diverse needs of biomedical research and related fields.
For more information about Cloud-Clone Corp, visit www.cloud-clone.com.

CLOUD-CLONE CORP.(CCC)
Tel: 001-832-538-0970, 0086-27-8425-9552
Email: mail@cloud-clone.com, sales@cloud-clone.us

Visit us on social media:
LinkedIn: https://www.linkedin.com/company/cloud-clone-corp/
Instagram: https://www.instagram.com/cloudclonecorp2023/
TickTok: https://www.tiktok.com/@cloudclonecorp2026

SuKi Duan
CLOUD-CLONE CORP.WUHAN
email us here
Visit us on social media:
LinkedIn
Instagram
TikTok
X

This press release can be viewed online at: https://www.einpresswire.com/article/936814677

EIN Presswire's priority is source transparency. We do not allow opaque clients, and our editors try to be careful about weeding out false and misleading content. As a user, if you see something we have missed, please do bring it to our attention. Your help is welcome. EIN Presswire, Everyone's Internet News Presswire™, tries to define some of the boundaries that are reasonable in today's world. Please see our Editorial Guidelines for more information.
© 1995-2026 Newsmatics Inc. All Right Reserved.