

To improve life freedom and happiness with Airwheel H3T electric automatic folding wheelchair.


To improve life freedom and happiness with Airwheel H3T electric automatic folding wheelchair.

BRUSSELS, BRUSSELS, BELGIUM, November 27, 2019 / EINPresswire.com/ -- Abstract: Airwheel H3T was born taking advantage of the trend and Airwheel will enable more users to experience the charm of technology and feel the freedom created by intelligence.

With the development of society and the improvement of human civilization, people, especially disabled people, are increasingly required to use modern high technology to improve their quality and freedom of life. While innovation has been an important part of Airwheel's legacy, our highest priority has always been customer satisfaction. So, Airwheel H3T automatic folding wheelchair was born and it is super easy to operate. The intelligent controller on the handlebar includes all basic functions, accelerating, decelerating, lighting, horns and folding, unfolding. The direction is controlled by the joystick, which ensures everyone can sit and go.



Airwheel H3T smart wheelchair



Airwheel H3T power wheelchair

The <u>Airwheel H3T electric wheelchair</u> is equipped with a Bluetooth remote

control, which remotely control the movement, folding and lighting of intelligent wheelchairs. It can help users to ride and store it by their own, with a higher degree of freedom. In addition, Airwheel also developed a mobile phone app, which can be used to control it through the mobile phone. It is the same control option as the Bluetooth remote control and smart controller.

The H3T medical equipment adopts an electric folding design. It only needs to control one button, and the body can be automatically unfolded and folded without laborious folding. It is especially suitable for people who have difficulty in walking, and can also bring convenience to the caregiver.

Safety is extremely important no matter for which product. The H3T power wheelchair uses the EABS electromagnetic brake. When the handle is released, it will brake immediately. Even on the slope, it can guarantee the parked state without slipping. When reversing, the horn will automatically sound a warning tone to remind surrounding passers-by and vehicles to avoid it in time. The rear side is also designed with a multifunctional anti-dumping wheel, which can effectively prevent the wheelchair from falling over when going up and down and backwards, and can also be used as an arm-brace when folding.



Airwheel H3T electric(power) wheelchair

Pls note the outstanding support and service we offer does not stop after delivery. You will always have free and unlimited access to our support team.

Airwheel Holding Limited Airwheel Holding Limited +32 10 61 33 24 email us here

This press release can be viewed online at: http://www.einpresswire.com

Disclaimer: If you have any questions regarding information in this press release please contact the company listed in the press release. Please do not contact EIN Presswire. We will be unable to assist you with your inquiry. EIN Presswire disclaims any content contained in these releases. © 1995-2019 IPD Group, Inc. All Right Reserved.