

Validity Diagnostics to Address Workplace Urine Drug Testing Subversion at SOFT-TIAFT 2026 Joint Meeting

Validity Diagnostics reagents are superior in detecting synthetic urine when it is submitted for a urine drug test

JACKSONVILLE, FL, UNITED STATES, August 4, 2026 /EINPresswire.com/ -- Validity Diagnostics (VDx) of Branford, FL, a leader in urine drug test validity solutions, will exhibit their products at the [SOFT-TIAFT 2026 Joint Meeting](#) in Chicago, IL on September 19-24, 2026. The company will highlight how advanced validity testing can help Fortune 500 companies, federal and state government employers, federal and state drug courts / probation offices and drug testing laboratories combat the growing threat of workplace urine drug test subversion.

As the submission of synthetic and adulterated urine continues to rise, many drug testing programs face a serious challenge - maintaining urine sample integrity. The Validity Diagnostics True Essentials Panel, which includes Oxidant History, True Urine LD and True Urine SD reagents, provides drug testing laboratories with dependable, science-backed tools to detect synthetic and/or adulterated urine.

"Workplace urine drug testing is only as dependable as the integrity of the urine sample" says CEO of Validity Diagnostics, Donna Suggs. "We are seeing increased sophistication in urine sample subversion and our products are specifically designed to help drug testing laboratories detect these threats before they compromise compliance or safety in the workplace."

Independent research, including the 2025 study by Franke, et. al. from Jena University Hospital, GERMANY and a 2025 study from [RTI International](#) in Research Triangel Park, North Carolina, demonstrated that the Validity Diagnostics True Urine LD and True Urine SD reagents outperformed conventional validity tests when identifying synthetic urine. These findings

APMC SYNTHETIC URINE CHALLENGE

9/23/2020

Sample	pH	Urine Creat	SGI	OXH	LD	SD	Result
Perfect Urine	7.50	11.8	1.0007	0.1	3	1.5	Fail
Dr Green Super Hero	7.30	69.6	1.013	1.6	9.3	3.4	Fail
Quick Fix Plus	7.30	67.3	1.0125	0.1	3	4	Fail
Ultimate Gold	6.40	39.9	1.0125	1.4	1.2	4.8	Fail
Ultra Pure	8.70	61.6	1.0117	0.5	8.4	0.4	Fail
Dr Greens X Agent	6.50	85.3	1.0102	2.3	9.3	4.3	Fail
Gator Whizz	6.80	77.1	1.0105	2.1	9.3	2.1	Fail
Golden Flask	5.50	43.8	1.014	3.3	9.3	5.9	Fail
Monkey Flask	5.00	44.7	1.0187	0.5	5.2	4.6	Fail
SyntheticS	8.10	35.1	1.0157	1.1	1.6	6.4	Fail
UPASS	6.70	78.8	1.0136	2.1	9.3	8.2	Fail

8/14/2025

Sample	pH	Urine Creat	SGI	OXH	LD	SD	Result
Test Clear PU100	7.99	98	1.0178	9.2	5.7	2.1	Fail
QuickFix Lot 1	6.96	45	1.0185	0	6.6	5.3	Fail
QuickLuck	6.00	31	1.0186	1.2	5.7	4.9	Fail
Pure Gold	5.60	7.8	1.0157	0	6.6	4.8	Fail
Xstream Lot #1	7.92	28	1.0161	0	3.9	4.1	Fail
Pure Steam	5.59	8.7	1.0154	0	4.8	5.6	Fail
Clean Whizz	9.38*	119	1.0125	0	3.9	5.4	Fail
Upass	9.88*	66	1.0144	0	4.8	4.6	Fail
Xstream Lot #2	7.01	32	1.0175	0	5.7	5.5	Fail
Monkey Whizz	5.59*	54	1.0161	1.8	5.6	2.4	Fail
Agent X	7.23	35	1.0194	0	2.1	4.1	Fail
QuickFix Lot 2	7.27	49	1.0185	0	3.9	6.9	Fail
Monkey Flask	9.45*	60	1.0049	0	4.8	-16	Fail
Golden Flask	7.63	19	1.0257	0	5.7	7.1	Fail
Sub Solution	7.56	81	1.0175	0	5.7	6.3	Fail

TAKE THE NEXT STEP AND CALL VDX TODAY

(855) 826-5500 | WWW.VDXCARES.COM

* Possible bacterial growth but no growth visible in tubes.
Validity Diagnostics SVT tests were performed on a Beckman Olympus AU400e in the toxicology laboratory at Advanced Pain Management Clinic, Jacksonville, FL.
Samples were received as unknowns labeled 1-11 in 2020 and 1-15 in 2025.




Synthetic Urine Detected by Validity Diagnostics reagents

confirm that traditional detection methods are no longer sufficient in today's urine testing environment where synthetic urine products continue to evolve.

About Validity Diagnostics

Validity Diagnostics (VDx) is an industry leader in laboratory innovation and validity testing for urine drug screening. The VDx True Essentials Panel, featuring True Urine LD, True Urine SD, and Oxidant History, helps laboratories and collection sites detect synthetic and adulterated urine samples with unmatched accuracy before the samples are tested for drugs of abuse or misuse. VDx partners with reference laboratories, hospital laboratories, and program administrators nationwide to promote urine drug testing integrity and reliability.



Learn more at www.vdxcare.com.



Workplace urine drug testing is only as dependable as the integrity of the urine sample."

Donna Suggs, CEO at Validity Diagnostics

Media Contact:

Ricky P. Bateh, PhD
Validity Diagnostics
rick@vdxcare.com
www.vdxcare.com

Ricky P Bateh, PhD
Validity Diagnostics
+1 904-448-3168
rick@vdxcare.com
Visit us on social media:
[LinkedIn](#)

This press release can be viewed online at: <https://www.einpresswire.com/article/856396633>

EIN Presswire's priority is source transparency. We do not allow opaque clients, and our editors try to be careful about weeding out false and misleading content. As a user, if you see something we have missed, please do bring it to our attention. Your help is welcome. EIN Presswire, Everyone's Internet News Presswire™, tries to define some of the boundaries that are reasonable in today's world. Please see our Editorial Guidelines for more information.

© 1995-2026 Newsmatics Inc. All Right Reserved.