

Cloud-Clone ELISA Kits: High-Precision Molecular Quantification Tools for Global Scientific Research

Complete ready-to-use ELISA systems with ultra-high sensitivity and cross-species compatibility to support preclinical and pharmaceutical research

HUSTON, TX, UNITED STATES, August 19, 2026 /EINPresswire.com/ -- In precision detection workflows across life science research, ELISA kits act as high-precision quantitative test strips to measure specific proteins, antibodies and small-molecule analytes within biological samples. With years of deep expertise in immunoassay development, [Cloud-Clone](#) has built proprietary core



technical platforms to manufacture premium ELISA kits covering thousands of targets across numerous species. Delivering accurate, consistent and user-friendly experimental solutions, these products are widely adopted in oncology, immunology, endocrinology, new drug development, food safety and a broad spectrum of research fields for scientists all over the world.

I. Cloud-Clone ELISA Kits: Fully Integrated, Ready-to-Use Complete Assay Systems

Designed under standardized all-in-one development principles, Cloud-Clone ELISA kits require no additional matching reagents and support immediate experimentation upon opening. Every kit contains a full set of essential components to form a closed, streamlined testing workflow:

1. Pre-coated microplate (core component): 96-well detachable plates pre-immobilized with high-specificity capture antibodies or antigens on well surfaces. The pre-coating eliminates manual coating errors from laboratory preparation, allows split testing across multiple sessions and cuts unnecessary reagent waste.

2. Detection antibodies and enzyme conjugates: Paired high-affinity detection antibodies combined with HRP and other enzyme labels form antigen-[antibody](#)-enzyme complexes that deliver robust signal amplification for colorimetric readouts.

3. Serial calibrator standards: 5 to 8 concentration gradients of purified target protein are supplied to establish broad-range standard curves, the fundamental basis for precise quantitative analysis.

4. Chromogenic and stop buffer system: TMB substrate solution (light-sensitive, turns blue under enzymatic catalysis) and stop solution (shifts the reaction to stable yellow color) deliver consistent color development with distinct concentration gradients.

5. Auxiliary consumables: Concentrated wash buffer to reduce matrix interference, plate sealing films to prevent cross-contamination, and detailed operation manuals. The full supporting package is compatible with serum, plasma, cell culture supernatant, tissue homogenate and other common biological matrices.

Figure 1 Cloud-Clone □Competitive Inhibition [ELISA Kit](#)

II. Four Core Assay Formats to Meet All Research Scenarios

Based on the lock-and-key principle of specific antigen-antibody binding, Cloud-Clone optimizes four mainstream ELISA methodologies tailored to distinct analyte properties and experimental objectives:

1. Sandwich ELISA (Most Widely Adopted Standard Format)

Adopting a capture antibody-target antigen-detection antibody sandwich structure, this format is optimized for large molecular antigens including cytokines, tumor biomarkers and immunoglobulins, exemplified by TNF α and BDNF detection kits. It features straightforward protocols completed within 3 hours, pg/mL-level ultra-sensitivity and exceptional specificity, with intra-assay CV below 10% and inter-assay CV below 15% to guarantee outstanding experimental reproducibility.

2. Competitive ELISA (Exclusive for Small Molecules)

Developed for low-molecular-weight substances such as hormones, neurotransmitters and drug residues that lack multiple antibody binding epitopes. The mechanism relies on competitive binding between sample analytes and labeled reference antigens, with signal intensity inversely correlated to sample concentration. With two decades of technical breakthroughs in small-molecule conjugation and immunogenicity tuning, Cloud-Clone's competitive ELISA portfolio fills critical gaps in global small-molecule quantification research, covering all hormone and metabolite targets.

3. Indirect ELISA (Dedicated to Antibody Detection)

Purified antigens are immobilized to capture target antibodies present in samples, followed by enzyme-conjugated secondary antibodies to amplify detection signals. Ideal for antibody screening and pathogen serology testing (such as HIV and viral antibody detection), this format delivers powerful signal amplification at controlled costs and suits high-throughput large-scale sample screening projects.

4. Direct ELISA (Rapid Preliminary Screening)

Streamlined workflows utilize enzyme-labeled antibodies to directly bind sample antigens, finishing all operations within just 2 hours for fast antigen qualitative analysis and preliminary antibody screening. Cloud-Clone's optimized formulation balances long-term reagent stability with baseline sensitivity to support quick preliminary compound and sample assessment.

III. Four Core Competitive Advantages Defining Premium-Grade ELISA Performance

1. Ultra-high sensitivity to detect trace biological signals

Through screening of high-affinity specific antibodies and iterative signal amplification technology, the limit of detection reaches the pg/mL range, enabling precise quantification of low-abundance analytes including trace hormones and cytokines. This resolves the common limitation of conventional kits that fail to measure dilute target molecules; for instance, the BDNF kit achieves sensitivity lower than 0.118 ng/mL, while the EPO kit accurately quantifies extremely low-concentration biological samples.

2. Superior specificity to eliminate cross-reactive matrix interference

Multiple rounds of stringent antibody screening and affinity refinement grant exceptional target recognition accuracy with negligible cross-reactivity against structural analogs including homologous hormones and paralogous proteins. Specialized matrix optimization for complex samples (serum, plasma, tissue homogenate etc.) maintains consistent recovery rates between 87% and 105%, effectively suppressing background interference.

3. Extensive cross-species coverage compatible with diverse sample types

The product library includes over 8,000 detection targets spanning dozens of species such as human, mouse, rat, bovine, porcine and rabbit to fit diverse animal model research demands. Kits support direct testing of serum, EDTA/heparin-anticoagulated plasma, cell supernatant, tissue homogenate, saliva and other matrices without pre-processing, drastically boosting laboratory operational efficiency.

4. Rigorous quality control for stable, publication-grade data

All manufacturing processes comply with ISO9001 and ISO13485 international quality management standards. Each production batch undergoes full verification of four core performance indicators: sensitivity, specificity, precision and long-term stability. Target activity attenuation remains under 5% throughout the shelf life, and standard curves deliver linear correlation coefficients $R^2 > 0.99$, generating data fully acceptable for SCI journal submissions. To date, more than 40,000 peer-reviewed research papers worldwide have cited Cloud-Clone ELISA kits.

IV. Four Key Selection Criteria to Optimize Kit Matching for Your Experiments

Researchers can avoid mismatched reagents and wasted trials by evaluating these four core indicators when selecting Cloud-Clone ELISA kits:

1. Confirm target analyte and species identity: Clarify the molecule type (protein, antibody or small molecule) and sample species. Cloud-Clone kits are strictly species-specific to prevent invalid cross-species testing.

2. Verify compatible sample matrices: Confirm whether the kit matches your sample type (serum, plasma, cell supernatant, etc.). Matrix-optimized kits are recommended for highly complex biological samples.

3. Review critical performance parameters:

- Sensitivity: Select pg/mL-grade kits for trace analytes, ng/mL-grade kits for routine target quantification.

- Detection range: Ensure the expected analyte concentration of your samples falls fully within

the kit's linear measurement interval.

- Precision: Prioritize products with intra-assay CV <10% and inter-assay CV <15%.

4. Opt for ready-to-use formulations: All Cloud-Clone ELISA kits come with pre-coated plates and liquid working buffers with no manual reagent preparation required, cutting human operation error and lowering the learning curve for novice laboratory staff.

Conclusion

From large intact proteins to trace small-molecule metabolites, from fundamental life science exploration to preclinical pharmaceutical evaluation, Cloud-Clone ELISA kits stand as trusted global detection tools built on ultra-sensitivity, outstanding specificity, consistent performance and broad compatibility. Moving forward, Cloud-Clone will continue advancing immunoassay technology, iteratively upgrading product performance and launching new kits tailored to cutting-edge research directions, driving breakthroughs and innovation across global life science communities.

SuKi Duan

CLOUD-CLONE CORP.WUHAN

[email us here](#)

Visit us on social media:

[LinkedIn](#)

[Instagram](#)

[TikTok](#)

[X](#)

This press release can be viewed online at: <https://www.einpresswire.com/article/935267042>

EIN Presswire's priority is source transparency. We do not allow opaque clients, and our editors try to be careful about weeding out false and misleading content. As a user, if you see something we have missed, please do bring it to our attention. Your help is welcome. EIN Presswire, Everyone's Internet News Presswire™, tries to define some of the boundaries that are reasonable in today's world. Please see our Editorial Guidelines for more information.

© 1995-2026 Newsmatics Inc. All Right Reserved.